

Airstel is a single-shot ready-to-use disinfectant designed for the bactericidal, mycobactericidal, fungicidal and virucidal disinfection of inaccessible surfaces. Its unique formulation incorporates a natural cinnamon oil extract which is especially effective against *Aspergillus brasiliensis*. Airstel is intended to be used for terminal disinfection, following contact surface disinfection, to ensure airborne microorganisms are eliminated.



APPLICATIONS

Airborne surface disinfection in:

- Operating theatres
- Intensive care units
- Other critical areas within the hospital
- Isolation wards
- Mortuaries

MICROBIOLOGICAL EFFICACY

Airstel is proven effective against microorganisms of concern such as:


- *Aspergillus brasiliensis*
- *Staphylococcus aureus*
- *Pseudomonas aeruginosa*
- Influenza A virus (H1N1)
- *Candida albicans*
- *Escherichia coli*
- *Mycobacterium avium*
- *Mycobacterium terrae*

HOW TO USE




Direct the valve away from you, press it down and turn it clockwise until it is locked.

Note: immediately after pressing down the lever, the Airstel formulation is released from the bottle.



Place the bottle on an even surface and leave the room. Airstel forms a cloud that covers approximately 40m³ in 1 minute (50ml bottle) and 150m³ in 3 minutes (300ml bottle).




Ventilate and allow a period of one hour before re-entering the premises.

Aspergillus mitospores

One of the challenges healthcare facilities face is with *Aspergillus* mitospores. They are relatively small (under 2.5 micron), which makes it difficult for air filters to capture them.

Aspergillus mitospores can settle in human lungs causing respiratory illnesses such as Aspergillosis. Effective airborne surface disinfection is essential.

Airstel's total discharge valve permits the nebulisation of the content in a single application. Within seconds a cloud of disinfectant is formed that will capture even the smallest of mitospores.



PRODUCT OPTIONS

Airstel 300ml

AG003



Airstel 50ml

AG001

Airstel is classified as a biocide under the EU Biocides Regulation. Use biocides safely. Always read the label and product information before use.

Airstel™