

TECHNICAL DATA SHEET

Ansell Encore Sterile Latex Gloves



DESCRIPTION

ENCORE® gloves are designed and engineered for surgeons, OR nurses and surgical technicians seeking protection and safety with advanced levels of ergonomics, dexterity and comfort during surgery. ENCORE® gloves afford the kind of hand protection that elevates users' confidence and provide optimum comfort during surgery. The Encore® Style 85 has been designed to meet the needs of users who want a glove with a rough surface for the sure handling of instruments, with maximum comfort and safety levels. Offering Ansell's SureFit™ technology to prevent cuff roll-down, this glove provides healthcare professionals with superior barrier protection.



TECHNICAL SPECIFICATION

Manufacture	Ansell
Shelf Life	3 Years
Size	S, M, L, XL, XXL
Quantity	Pack of 100 (50 pairs)
Material	Latex with polyurethane inner coating
Test Standards	EN 455 EN 374 (parts 1, 2, 4 and 5) EN 420 EN 556 ASTM D3577 AST D6978 ASTM D7160

FEATURES

- Textured, Chlorinated Finish
- Ergonomically Designed for Comfort
- Advanced Dexterity
- PU (Polyurethane) Inner Coating Resistant to Alcohol-Based Disinfectants
- Low Protein Content for Minimal Skin Irritation
- Natural Latex
- Powder-Free