

TECHNICAL DATA SHEET



Berkshire Isolator Tool Padded Mop Covers

DESCRIPTION

Efficient cleaning of compounding isolators, biological safety cabinets, laminar flow hoods and other mini environments that are critical to the quality and purity of compounded products. The padded knitted cover and rounded edges of the mop head ensure the tool does not scratch or damage surfaces.

Mop covers can be Gamma Irradiated but are supplied non-sterile. Covers are autoclavable.



TECHNICAL SPECIFICATION

Quantity	Pack of 25 covers
Compatibility	ISO class 4-10 cleanroom
Sterility	Non-sterile

FEATURES

- 100% continuous filament knitted padded covers are laundered and packaged in an ISO Class 4 (Class 10) cleanroom
- Foam-free cushioned design insures even surface coverage with reduced residues and contaminants
- Easy-to-install covers fit securely and stay in place during use
- Autoclavable