

TECHNICAL DATA SHEET



BioClean Chemotherapy Preparation Mat Non-Sterile

DESCRIPTION

The BioClean ChemPrep™ Prep Mat enhances workbench safety for technicians, nurses and pharmacists when handling or administering chemotherapy drugs. The Prep Mat is a highly-absorbent, three-layered preparation mat consisting of two layers of polyethylene with a central layer of air-laid paper. Suitable for use in ISO Class 5 or EU GMP Grade A cleanroom. Each mat is double bagged.



TECHNICAL SPECIFICATION

Size	Each mat is 400mm (16") x 250mm (10")
Weight	250gsm
Compatibility	ISO Class 5 or EU GMP Grade A cleanroom
Quantity	Pack of 10
Colour	Blue

FEATURES

- Protects against spillages
- Highly-absorbent
- Latex-Free 3 layer construction
- Suitable for ISO Class 5 cleanrooms
- 250gsm weight
- Low-linting
- Double-bagged

Each mat is made up of three layers:

- Top Layer - Polyolefin apertured film evenly distributes liquid
- Middle Layer - Highly absorbent air-laid paper
- Bottom Layer - Polyethylene defends against leakages