

TECHNICAL DATA SHEET

Biohazard Spill Kit



DESCRIPTION

Multiple use kit with tools and PPE for safe clear up and disinfection of multiple biohazard spillages. NaDCC chlorine release granules absorb spillage, kill pathogens on contact and reduce cross contamination risk. The tablets dilute to produce 10,000ppm concentration chlorine disinfectant solution fully effective against blood borne pathogens and Tuberculosis (TB). Intended for use in busy emergency rooms, larger departments and other similar situations where the risk and frequency of biohazard spillages is high. Refill items can be purchased separately.



TECHNICAL SPECIFICATION

Quantity	Mini bodily fluids kit (2 spills) Standard bodily fluids kit (4 spills) Midi biohazard kit (10 spills) Multi biohazard kit (25 spills)
Shelf Life	3 years (some kit components may vary)
Usage	NHS hospitals, clinics, surgeries, laboratories, ambulances, dental and veterinary practices, nursing and care homes, schools, offices, businesses, ships and aircrafts
Test Standards	EN 1650 EN 1276 EN 13624 EN13727 EN 13697

FEATURES

- Portable and Lightweight
- Multiple uses form a single kit
- Safe
- Efficient
- Conforms to WHO and CDC guidelines
- Ready to use
- Easy to keep
- PPE included
- Refills available

TECHNICAL DATA SHEET

Biohazard Spill Kit– Mini (2 spills)



Mini Kit Contents:

1 x Blood/ Biohazard pack

1 x Urine/ Vomit pack

The standard bodily fluid kit contains 2x of each pack. The midi and multi biohazard kits contain enough supplies to clean 10 and 25 spills respectively.

TECHNICAL DATA SHEET

Biohazard Spill Kit– Standard (4 Spills)



Blood/Biohazard Pack Contents:

- 1 x pair nitrile gloves
- 1 x apron
- 1 x 100g SoChlor NaDCC / tablet diluter
- 1 x clinical waste bags
- 4 x SoChlor TAB chlorine tablets
- Paper towels
- Illustrated instructions on pack

Urine/Vomit Pack Contents:

- 1 x pair of nitrile gloves
- 1 x apron
- 1 x sachet of GV High Absorbency Granules
- 1 x scoop and scraper
- 2 x packs of disinfectant surface wipes
- 1 x Clinical waste bag
- Illustrated instructions on pack

TECHNICAL DATA SHEET

Biohazard Spill Kit– Midi (10 Spills)



Kit Contents:

- 100 x 2.5g Chlorine tablets
- 1 x 500g NaDCC Granules
- 1 litre diluter
- 10 pairs nitrile gloves
- 10 aprons
- 10 clinical waste bags
- 10 x scoops and scrapers
- Illustrated instructions / poster
- Recording card
- Hard sided, lightweight bright yellow case
- Wall mount

TECHNICAL DATA SHEET

Biohazard Spill Kit– Multi (25 Spills)



Kit Contents:

- 100 x 2.5g Chlorine tablets
- 1 x 500g NaDCC Granules
- 1 litre diluter
- 25 pairs nitrile gloves
- 25 aprons
- 25 clinical waste bags
- 25 x scoops and scrapers
- Illustrated instructions / poster
- Recording card
- Hard sided, lightweight bright yellow case
- Wall mount

TECHNICAL DATA SHEET

Biohazard Spill Kit



Approved for use on the following

- Blood
- Breast Milk
- Amniotic Fluid
- Vaginal Secretions
- Semen
- Cerebrospinal Fluid (CSF)
- Synovial Fluid
- Pericardial Fluid
- Peritoneal Fluid
- Unfixed Tissues/Organs
- Saliva