

TECHNICAL DATA SHEET

Byotrol Trionic Biocide Spray



DESCRIPTION

A non-alcohol formulation for the decontamination of medical devices and equipment, cleaning surfaces (including the removal of bio-film), and then disinfecting with broad spectrum biocidal efficacy. The disinfectant kills bacteria, fungi, yeasts, viruses and mycobacteria, and has been tested to be effective under a number of microbiological test conditions.

Note: This product has been discontinued by the manufacturer. We do have a small amount of stock available, so please contact us to check quantities before placing an order.



TECHNICAL SPECIFICATION

Quantity	750ml Trigger Spray
Test Standards	EN 1040 (Bactericidal) EN 1276 (Bactericidal) EN 13727 (Bactericidal) EN 1650 (Fungicidal) EN 13697 (Bactericidal, Yeasticidal, Fungicidal) EN 13624 (Fungicidal, Yeasticidal) EN 14476 (Fungicidal) EN 14348 (Mycobacterial)
Recommended Use	Medical devices in clinical and other human health environments, as well as cleanrooms and laboratories.

FEATURES

- Alcohol-Free
- Quat-Based Biocide
- Cleans and Disinfects
- Ready to Use (RTU)
- Proven effective against a wide range of micro-organisms
- Safe and convenient to use