

TECHNICAL DATA SHEET



Chemgene High Level Disinfectant Wipes

DESCRIPTION

High level quaternary ammonium based disinfectant wipes. Manufactured from a nonwoven material, offering strength, absorbency and low lint properties at a cost effective price. Wipes are impregnated with a multi-active disinfectant solution that is alcohol free. Recommended for the cleansing and disinfecting of hard surfaces of non-invasive, non-lumened medical devices such as monitors, IV poles and X-Ray Equipment.

Please note: Tubs of 200 are CE marked as a class IIa medical device.



TECHNICAL SPECIFICATION

Sterility	Non-sterile
Quantity	100 wipes per pouch 200 wipes per tub
Test Standards	EN 13727 EN 13624 EN 13697 EN 16615 EN 14476

FEATURES

- Powerful disinfectant
- Non woven wipes
- Alcohol free
- Effective against a range of micro-organisms
- Non-corrosive

TECHNICAL DATA SHEET



Chemgene High Level Disinfectant Wipes

This product has been tested effective against:

Virucidal Activity

- Adenovirus
- Herpes Simplex
- Vaccina Virus
- Hepatitis B
- Murine Norovirus
- Enveloped viruses (including SARS CoV-2)

Bactericidal Activity

- Candida albicans
- Escherichia coli
- Methicillin-resistant Staphylococcus aureus
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Enterococcus hirae
- Listeria monocytogenes
- Methicillin-resistant Staphylococcus pseudomonas
- Salmonella tryphimurium