

TECHNICAL DATA SHEET



Chemgene Instrument Disinfectant Wipes

DESCRIPTION

High level non-woven disinfectant wipes for the cleaning of non-surgically invasive instruments, such as transvaginal and transrectal probes. Also suitable for use on hard surfaces of non-invasive, non-lumened medical devices, such as monitors, radiographic and diagnostic equipment. The wipes are impregnated with a multi-active disinfectant solution that is alcohol free and tested to EN standards as effective against a range of bacteria, viruses, fungi, yeasts and mycobacteria. CE marked as a class IIb medical device. Not to be used to disinfect surgically invasive medical devices.



TECHNICAL SPECIFICATION

Quantity	100 wipes per tub
Test Standards	EN 14776 EN 1650 EN 13697 EN 1276

FEATURES

- Hard surfaces
- Non-surgically invasive instruments
- Non-invasive, non-lumened medical devices
- Powerful disinfectant
- Non-woven wipes
- Alcohol free
- Non-corrosive
- Class IIb CE marked
- Effective against a range of micro-organisms

TECHNICAL DATA SHEET



Chemgene Instrument Disinfectant Wipes

This product has been tested effective against:

Virucidal Activity

- Hepatitis B
- Herpes Simplex
- Norovirus Vaccinia
- Vaccinia
- Adenovirus

Bactericidal Activity

- Pseudomonas aeruginosa
- Staphylococcus aureus
- Legionella pneumophila
- Listeria monocytogenes
- Methicillin-resistant Staphylococcus aureus
- Escherichia Coli
- Enterococcus hirae
- Salmonella typhimurium
- Methicillin-resistant Staphylococcus pseudointermedius

Bacteriophage Activity

- Lactococcus lacti

Fungicidal Activity

- Candida albicans
- Trichophyton rubrum