

TECHNICAL DATA SHEET



Chemgene MedLab Non-Sterile Surface Disinfectant

DESCRIPTION

A high level disinfectant for surfaces in a range of applications, including laboratory, forensics and healthcare, to prevent the transmission of harmful pathogens. Chemgene MedLab Multi-Surface disinfectants are suitable for the cleaning and disinfection of walls, floors, surfaces, instruments and equipment. Supplied in 1L and 5L concentrate formats, with the options of eucalyptus fragrance or unfragranced.

Note: This product replaces the Chemgene Medical and Chemgene Laboratory ranges of disinfectants.



TECHNICAL SPECIFICATION

Quantity	6 x 1L bottles 2 x 5L bottles
Scent Options	Unscented Eucalyptus
Test Standard	EN 13727 (bactericidal) EN 17387 (bactericidal) EN 13624 (yeasticidal, fungicidal) EN 17387 (yeasticidal, fungicidal) EN 16777 (virucidal - all enveloped viruses) EN 16777 (Adenovirus, Norovirus) EN 14476 (virucidal - all enveloped viruses) EN 14476 (Adenovirus, Norovirus) EN 13623 (Legionella) EN 13610 (bacteriophage)

FEATURES

- Non-Corrosive
- Quicker than comparable disinfectants
- Suitable for use on bio-hazard contamination
- Tested against a range of standards
- Available in 1L and 5L formats