

TECHNICAL DATA SHEET

Chemo Preparation Mat— Non-Sterile



DESCRIPTION

The disposable non-sterile chemotherapy prep mat provides an ideal surface for working with chemotherapy and cytotoxic drugs. The surface of the 3-layer mat is firm and flat with a penetrable mesh top layer. The second layer is absorbent gauze, which absorbs spillages, and the bottom layer is non-permeable plastic, preventing spills from leaking through to work surfaces.



TECHNICAL SPECIFICATION

Quantity	Case of 50 mats
Sterility	Non –Sterile
Dimensions	559mm x 406mm (22" x 16")

FEATURES

- Disposable
- Non-sterile
- 3 layers
- Non-linting
- Single Use