

TECHNICAL DATA SHEET



Chemotherapy Protective Apron - Sleeveless

DESCRIPTION

The Chemotherapy Protective Apron with Sleeves provides protection against a range of cytotoxic drugs. (Full test data sheet and certification available on request) . Manufactured from lightweight low linting material, the apron features an adjustable neck fastening, tie tapes at the rear and ultrasonically bonded sleeve seams with protective tape. Tested against permeation standard ASTM F739-12. Ultrasonically bonded seams with protective tape. Silicone free.



TECHNICAL SPECIFICATION

Material	CleanTough
Construction	Adjustable neck, tie fastening at waist
Classification	PPE Cat 3 type PB (6)
Compatibility	ISO Class 4
Characteristics	Low Particulating
Sterility	Non-sterile (Sterile available on request)
Size	Medium Large Extra Large

FEATURES

- Ultrasonically bonded sleeve seams with protective tape & 100% polyester elasticated cuffs
- PPE Cat 3 Type PB [6]
- Low participating
- ISO Class 4
- Also available sterile
- For single use
- Lightweight and hard wearing