

TECHNICAL DATA SHEET

Chemotherapy Protective Apron



DESCRIPTION

The Chemotherapy Protective Apron with Sleeves provides protection against a range of cytotoxic drugs. (Full test data sheet and certification available on request) . Manufactured from lightweight low linting material, the apron features an adjustable neck fastening, tie tapes at the rear and ultrasonically bonded sleeve seams with protective tape. Tested against permeation standard ASTM F739-12. Ultrasonically bonded seams with protective tape. Silicone free.



TECHNICAL SPECIFICATION

| | |
|-----------------|---|
| Material | Polypropylene, Polyethylene |
| Construction | Adjustable neck, tie fastening at waist |
| Size | Medium Large Extra Large |
| Characteristics | Low Particulating |
| Compatibility | ISO Class 4 |
| Classification | PPE Cat 3 type PB (6) |

FEATURES

- Ultrasonically bonded sleeve seams with protective tape & 100% polyester elasticated cuffs
- PPE Cat 3 Type PB [6]
- Low Particulating
- ISO Class 4
- Also available sterile
- For single use
- Lightweight