

# TECHNICAL DATA SHEET

## Chemotherapy Protective Oversleeves



### DESCRIPTION

The Chemotherapy Protective oversleeves provide additional protection against a range of cytotoxic drugs. (Full test data sheet and certification available on request). Manufactured from lightweight low linting material, the sleeve covers are UNIVERSAL SIZE. Tested against permeation standard ASTM F739-12. Ultrasonically bonded seams with protective tape. Silicone free. For Single Use. Also available sterile - please enquire.



### TECHNICAL SPECIFICATION

Construction	Ultrasonically bonded sleeve seams with protective tape & 100% polyester elasticated
Quantity	Pack of 30 sleeves
Classification	PPE Cat 3 Type PB [6]
Compatibility	ISO Class 4
Characteristics	Low participating

### FEATURES

- Low lint material
- Universal size
- Ultrasonically bonded seams with protective tape
- Silicone free
- 100% Polyester elasticated cuffs
- Additional protection
- For single use