

TECHNICAL DATA SHEET

Cleanroom Wipes Non-Sterile

100% Polyester



DESCRIPTION

Cleanroom Supplies Non-Sterile wipes are a superior quality, ultra low lint clean room dry wipe. This laser edge sealed Polyester wiper has high absorbent properties combined with cleanliness and strength. Suitable for low-abrasive applications.



TECHNICAL SPECIFICATION

Manufacture	Optimum Protection
Material	100% single Knit Polyester
Weight	135 g/m ²
Wipe Size	23cm x23cm (9" x 9")
Quantity	Pack of 50 or carton of 1500
Compatibility	GMP Grade A areas and ISO class 4-8 cleanrooms (FED 10 to 100,000)
Sterility	Non-Sterile

FEATURES

- Absorbs and retains solvents and chemicals
- Superb wet strength
- Fully laser sealed border to avoid particle release
- High abrasive and dissolvability resistance
- Ultra Low Lint
- Soft Texture
- Double Bagged
- No detectable silicone, amide or DOP