

TECHNICAL DATA SHEET

Cytotoxic Spill Kit



DESCRIPTION

Multi-use, portable and lightweight spill kit for intensive and frequent use on multiple and/or large scale spillages of cytotoxic drugs. High absorbency granules absorb and contain liquid spills. Absorbent pads and sterile water mop up powder spills. A custom solution is provided to neutralise cytotoxic residue of spill site and the surrounding area. Warn sign and colour coded, sealable bags are provided to aid safe disposal. Supplied gloves comply with ASTM D6978-05: Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.



TECHNICAL SPECIFICATION

Recommended Use	Used throughout the NHS in hospitals, clinics, surgeries, out-patient departments, pharmacies, laboratories; nursing and care homes and other locations where cytotoxic or chemotherapy drugs are handled and stored. Safe for use by non-specialist staff.
Test Standard	ASTM D6978-05 (Protective Gloves)

FEATURES

- Portable and lightweight
- Safe
- Efficient
- Conforms to WHO guidelines
- Ready to use
- Easy to keep
- Single use packs
- PPE included

TECHNICAL DATA SHEET



Cytotoxic Spill Kit

PPE Pack (x2) Contents:

- 3x Pairs of disposable gloves*
- 3x Pairs of oversleeves
- 3x Thumb loop gowns
- 3x Pairs of overshoes
- 3x FFP3 Masks
- 1x Kneeling pad
- 3x Eye visors

Gloves comply with ASTM D6978-05: Standard practice of assessment of resistance of medical gloves to permeation by chemotherapy drugs

Liquid Clean-Up:

- 3x Cyto granule 50g sachet
- 3x Scoop scraper
- 1x Pair of tweezers

Powder Clean Up:

- 3x Sterile water pods
- 3x Absorbent pads

Decontamination:

- 3x Cytali powder 10g sachet
- 1x 250ml Diluter
- Paper towels

Waste Disposal:

- 1x Instruction sheet
- 3x Cytotoxic waste bags
- 3x Bag ties
- 1x Warning sign
- 1x Record card