

TECHNICAL DATA SHEET

Disposable Long Sleeve Plastic Gown



DESCRIPTION

A disposable gown with an apron-style open back, long sleeves and thumb loop cuffs. Manufactured from an impervious PE fluid resistant film, making it suitable for use against prolonged contact with fluids. Suitable for use in medical industries by both staff members and patients.



TECHNICAL SPECIFICATION

Sterility	Non-Sterile
Colour	Blue
Quantity	Pack of 15 aprons Carton of 10 x flat aprons

FEATURES

- 43" (110cm) in Length
- 25 Micron Thickness
- Fluid Resistant
- Thumb Loops and Long Sleeves
- Latex Free
- Category I PPE
- Universal Size
- Single Use