

TECHNICAL DATA SHEET

Disposable Surgical Isolation Gowns



DESCRIPTION

Single-use, SMS isolation gown for protection against bodily fluids, micro-organisms and contaminated solids. 100% latex-free and suitable for surgical procedures with low-to-medium fluid exposure. The 40gsm SMS material is liquid repellent, breathable and low-linting. Cuffs are manufactured from knitted polyester for maximum comfort. A velcro neck fastening further increases protection and the waist belt can be tied around the back to provide an improved fit.



TECHNICAL SPECIFICATION

| | |
|---------------|---------------------------------------|
| Material | 100% Polypropylene |
| Quantity | Pak of 10 or carton of 100 |
| Sterility | Non-sterile |
| Shelf Life | 3 years from manufacture |
| Test Standard | EN 13795 (Surgical Clothing & Drapes) |

FEATURES

- 40gsm SMS Material
- CE Certified as a Class I Medical Device
- Single-Use
- Knitted Polyester Cuffs
- Low-to-Medium Level Fluid Protection
- Latex-Free
- Liquid repellent, breathable and low linting
- Velcro Neck Fastening
- Waist Belt for Improved Fit