

TECHNICAL DATA SHEET

Disposable Surgical Isolation Gowns



DESCRIPTION

A single-use sterile isolation gown which offers an affordable solution for procedures where low to moderate levels of fluid are expected. Constructed from 48gsm SMMS fabric, it provides high level barrier properties while also offering a high level of breathability to ensure optimal comfort for the wearer. The ergonomic fit with raglan sleeves provides the wearer with added comfort, a greater range of movement and allows for easier donning. A velcro neck fastening further increases protection and the waist ties (2x inner and 2x outer) provide an improved fit.



TECHNICAL SPECIFICATION

Material	100% Polypropylene
Quantity	Carton of 40 gowns
Size	M L XL XXL
Sterility	Sterile
Shelf Life	3 years from manufacture
Test Standard	EN 13795 (Surgical Clothing & Drapes)

FEATURES

- 48gsm SMMS Material
- CE Certified as a Class I Medical Device
- Single-Use
- Individually Wrapped
- Knitted Polyester Cuffs
- Low-to-Medium Level Fluid Protection
- Latex-Free
- Liquid repellent, breathable and low linting
- Velcro Neck Fastening
- Multiple Waist Ties for Improved Fit