

TECHNICAL DATA SHEET

Distel Disinfectant Dosing System



DESCRIPTION

A concentrated disinfection dosing system intended for the disinfection of hard, non-porous medical surfaces within the hospital environment. Every dose is the same volume as the previous one. It destroys DNA and RNA. The bottle heads of the Distel Dosing System make certain that each fresh dose of surface disinfectant is the same volume. Simply invert the dosing bottle once to prime and a second time to apply 5ml of surface disinfectant. Each dose is consistent, resulting in uniform disinfection for increased compliance.



TECHNICAL SPECIFICATION

Quantity	1 x 1L concentrate disinfectant, 1 x 1 Litre dosing bottle also available
Recommended Use	Monitors, keyboards, hoists, handrails, patient trolleys and other non-invasive medical devices and hard surfaces

FEATURES

- Reusable bottles
- 1 x 1L bottles of concentrate disinfectant produces 2000 doses
- Can be used with any dry wipe
- Fragrance Free
- Sold in 1L bottles only
- The Distel Dosing System is CE Marked as a Class IIa Medical Device.
- Effective against bacteria, yeast and viruses (including coronavirus) from 3 minutes