

TECHNICAL DATA SHEET

Infrared Forehead Thermometer



DESCRIPTION

A HIGH QUALITY and accurate infrared thermometer which measures the temperature of the human body or an object based on the infrared energy emitted by the object. As well as body temperature, this thermometer can be used to measure the temperature of liquids such as food or drink. Simply point the temperature probe at the target to get a quick measurement. In forehead mode, the measuring range is between 32.0°C - 42.9°C (89.6°F - 109.2°F). In object mode, measuring range is 0°C - 100°C (32°F - 212°F). The thermometer is able to store the last 35 readings, allowing you to track changes in temperature.



TECHNICAL SPECIFICATION

Recommended Use	Forehead, food and fluids
Warranty	2 years from manufacture
Quantity	1 thermometer, contact us if you would like to bulk order

FEATURES

- Easily switch between °C and °F
- No-contact infrared thermometer
- 2 modes - forehead and object temperature
- Auditory and visual fever warning
- Sounds can be switched on and off
- Memory stores 35 readings
- Forehead/body measuring range between 32.0°C - 42.9°C (89.6°F - 109.2°F)
- Object measuring range between 0°C - 100°C (32°F - 212°F)
- Supplied with 2 x AAA batteries and a storage bag