

TECHNICAL DATA SHEET

Micronclean Delta Sterile Disinfectant



DESCRIPTION

Micronclean Delta sterile neutral detergent, a hard surface cleaning agent designed for use in cleanrooms prior to disinfection. It is pH neutral after irradiation, helping to protect delicate cleanroom surfaces and those prone to corrosion. Manufactured with Water for Injection (WFI) quality water, for high purity. Safe to use on any surface. This product is filled and capped in ISO Class 5 uni-directional airflow cabinets, and packed within an ISO class 7 cleanroom. The disinfectant is filtered through a 0.2 micron filter at the point of filling and the product is gamma irradiated with a minimum dose of 25kGy.



TECHNICAL SPECIFICATION

Sterility	Sterile
Format	RTU Trigger Spray 15 x 950ml RTU 2 x 5L bottle
Filtration Level	0.2 micron filtered
Compatibility	ISO Class 7 cleanroom
Irradiation Standard	Gamma (25kGy)

FEATURES

- pH Neutral Detergent
- Sterile (Gamma Irradiated)
- Colourless
- Odourless
- All Formats Ready-To-Use (RTU)
- For Rotation with Micronclean Alpha and Beta