

TECHNICAL DATA SHEET

Non Sterile 70% DE



DESCRIPTION

70% Denatured ethanol alcohol (DE) in water for injection quality water (WFI). All products use European Pharmacopeia grade alcohol blended with WFI quality water to guarantee endotoxin levels below 0.25 EU/ml. Our 70/30 DE alcohol has many uses and is manufactured to the highest standard. In the cleanroom setting, alcohol disinfectants refer to two chemical compounds, denatured ethyl alcohol (DE) and isopropyl alcohol (IPA) also referred to as isopropanol. IPA and DE are broad spectrum antimicrobial and are rapidly bactericidal, fungicidal and viricidal, killing most bacteria within five minutes of exposure.



TECHNICAL SPECIFICATION

Irradiation Standard	Non-Sterile
Product	0.2 micron filtered
Manufacture	Agma
Sizes Available	1L and 5L bottles
Usage	Pharmaceutical and hospital cleanrooms
Chemical Formulation	70% IPA with water

FEATURES

- Fully assembled and ready for use
- 2 years unopened shelf life from manufacture
- 0.2 micron filtered
- Low endotoxin level– suitable for product contact areas
- Available in 900ml T/S and high quality wet wipes
- Gamma irradiated for sterility ~(20-40kGy)
- WFI Quality Water– High Purity Water
- Double Bagged