

TECHNICAL DATA SHEET

Non-Sterile Neutral Detergent



DESCRIPTION

A non-sterile, neutral detergent specifically developed for cleanroom use prior to disinfection. With highly effective surface cleaning properties, it is ideal for deep cleaning prior to disinfection, to remove visible soiling or build up of residue. Formulated to be low foaming, it has low residue and smearing properties. Blended with water for injection quality water (WFI), manufactured during a process guaranteeing a specification of less than 10 cfu per 100ml of water.

All products are available with full Quality Assurance documentation.



TECHNICAL SPECIFICATION

Irradiation Standard	Non-sterile
Product	0.2 micron filtered
Sizes	1L and 5L bottles
Quantity	15 x 1L bottles 4 x 5L bottles
Compatibility	Grade A conditions in a Grade C (ISO 7) cleanroom
Shelf Life	2 years from manufacture

FEATURES

- Ideal for dosing dry wipes and mop covers
- Bottles are individually double bagged
- 2 different bottle size options are available for small and large areas
- Low Endotoxin level – suitable for product contact areas
- Filled under Grade A conditions in a Grade C (ISO 7) cleanroom