

TECHNICAL DATA SHEET

CPE/PP Plastic Anti-Slip Shoe Cover



DESCRIPTION

Heavy duty CPE sole with good and effective anti-slip properties due to micro-indentations. Elasticated with polypropylene upper and CPE coated sole. Water resistant and heavy duty. Thicker gauge material tank standard CPE overshoes.



TECHNICAL SPECIFICATION

Manufacture	Optimum Protection
Material	Polypropylene / CPE
Sizes	Regular– size 5-10 shoes Large– size 10-14 boots
Quantity	Regular– 500 shoes (250 pairs) per carton Large– 500 shoes (250 pairs) per carton

FEATURES

- Water resistant sole
- Thicker Gauge Material than standard CPE Overshoes
- Elasticated
- Heavy duty
- With polypropylene upper and CPE coated sole