

TECHNICAL DATA SHEET



STAINLESS STEEL PETRI DISH STAND

Single-Tier Circular

DESCRIPTION

A high-quality, self-assembly stainless steel floor stand for cleanroom settle plates/petri dishes. EU GMP Annex 1 requires Cleanroom Environmental Monitoring to analyse the efficiency of your cleanroom disinfection and rotational disinfection regime. This petri dish stand is ideal for positioning plates in a non-intrusive but optimal location for cleanroom environmental monitoring. Mounted on a single rod with a solid, circular base manufactured from 6mm thick stainless steel for improved stability. The top contains 4 upturned tabs to prevent dishes from sliding off if the stand is disturbed.

Designed to your exact specifications.



TECHNICAL SPECIFICATION

Product	Stainless Steel Petri Dish Stand
Manufacturer	Optimum Protection™
Standard Material	304 Stainless Steel
Standard Finish	DP1 (Brushed)
Other Options	316 Stainless Steel Adjustable Height (Telescopic) Electropolishing Custom Sizes

FEATURES & USES

- 304 Grade Stainless Steel (316 Also Available)
- DP1 Finish As Standard
- Bespoke Design to Your Requirements
- Electropolishing Available
- Solid Baseplate for Improved Stability
- Designed to Hold 1x Plate/Dish
- Upturned Lips Prevent Dishes From Falling Off
- Allows for Non-Intrusive Environmental Monitoring

ADDITIONAL IMAGES



Height Adjustable
Stand



Stable Base Plate



Electropolished Finish