

# TECHNICAL DATA SHEET



## STAINLESS STEEL PETRI DISH STAND

### 3-Tier

#### DESCRIPTION

A high-quality, self-assembly stainless steel floor stand for cleanroom settle plates/petri dishes.

EU GMP Annex 1 requires Cleanroom Environmental Monitoring to analyse the efficiency of your cleanroom disinfection and rotational disinfection regime. This petri dish stand is ideal for positioning plates in a non-intrusive but optimal location for cleanroom environmental monitoring. Mounted on a single rod with a solid, circular base manufactured from 6mm thick stainless steel for improved stability.

Designed to your exact specifications.



#### TECHNICAL SPECIFICATION

<b>Product</b>	Stainless Steel Petri Dish Stand
<b>Manufacturer</b>	Optimum Protection™
<b>Standard Material</b>	304 Stainless Steel
<b>Standard Finish</b>	DP1 (Brushed)
<b>Other Options</b>	316 Stainless Steel Adjustable Height (Telescopic) Electropolishing Custom Sizes

#### FEATURES & USES

- 304 Grade Stainless Steel (316 Also Available)
- DP1 Finish As Standard
- Bespoke Design to Your Requirements
- Electropolishing Available
- Solid Baseplate for Improved Stability
- Designed to Hold 1x Plate/Dish and 1x Lid
- Upturned Lip Prevents Dishes From Falling Off
- Allows for Non-Intrusive Environmental Monitoring

#### ADDITIONAL IMAGES



Height Adjustable  
Stand



Stable Base Plate



Electropolished Finish