

TECHNICAL DATA SHEET

Stainless Steel Straight Sided Jug



DESCRIPTION

Stainless steel straight sided jug, manufactured from quality grade 304 stainless steel, this product can also be manufactured in grade 316 stainless steel. It comes fully welded and crevice free and is suitable for use in cleanrooms, laboratories and other clean environments. Bright polished and various options for different capacities.



TECHNICAL SPECIFICATION

Manufacture	Optimum Protection
Capacity	0.5L 1L 1.5L 2L

FEATURES

- 304 Grade stainless steel
- Can be manufactured in 316 grade stainless steel
- Fully welded
- Crevice free
- Bright polished finish
- Various capacities available