

TECHNICAL DATA SHEET

Stainless Steel Table & Undershelf



DESCRIPTION

Fully welded stainless steel construction for cleanroom, laboratory and food preparation areas. These products are designed and manufactured with quality and strength in mind. Ideal for sterile and hygienic environments, they provide a hygienic and easily cleanable working surface. The legs are manufactured from stainless steel box section, and have low level box section tie bars on 3 sides, allowing under bench units to be placed underneath. Each leg is fitted with adjustable feet which are used to even up the table level on un-even floor surfaces.



TECHNICAL SPECIFICATION

Steel Grade	304 or 316
Rear Upstand	Optional
Sizes	1200mm L x 600mm D x 900mm H 1500mm L x 600mm D x 900mm H 1800mm L x 600mm D x 900mm H 2100mm x L 600mm D x 900mm H
Recommended Use	Cleanrooms, laboratories, aerospace, defence, electronics, food processing, pharmaceutical processing, biotechnological processing and computer assembly

FEATURES

- Fully welded
- Hygienic
- Strong
- High quality
- Easy to clean
- 304 or 316 grade stainless steel