

# TECHNICAL DATA SHEET



## Stainless Steel Table With Castors – Cleanroom Supplies

### DESCRIPTION

Fully welded stainless steel construction for cleanroom, laboratory and food preparation areas. These products are designed and manufactured with quality and strength in mind. Ideal for sterile and hygienic environments, they provide a hygienic and easily cleanable working surface. Suitable for hygienic use within cleanrooms and laboratory environments. Tabletops are manufactured from 304 Grade stainless steel, with Grade 430 legs. All welds are ground and polished to ensure minimum contamination pockets. Our furniture is also available in higher grade 316 grade stainless steel.



### TECHNICAL SPECIFICATION

Stainless Steel Grade	304 with grade 430 legs
Sizes	1200mm L x 600mm D x 900mm H 1500mm L x 600mm D x 900mm H 1800mm L x 600mm D x 900mm H 2100mm x L 600mm D x 900mm H
Recommended Use	For cleanroom, laboratory, food preparation areas and pharmaceutical

### FEATURES

- Fully welded
- High quality
- Strong
- Made from grade 304 stainless steel
- Can be manufactured in grade 316 stainless steel