

TECHNICAL DATA SHEET



Stainless Steel Table With Castors – Electropolished Top

DESCRIPTION

Fully welded stainless steel construction for cleanroom, laboratory and food preparation areas. These products are designed and manufactured with quality and strength in mind. Ideal for sterile and hygienic environments, they provide an easily cleanable working surface. Suitable for hygienic use within cleanrooms and laboratory environments. Electropolished surfaces greatly reduce bacterial colonisation. The legs are manufactured from fully welded 30mm stainless steel box section, and have low level box section tie bars. Each leg is fitted with lockable castors to enable easy transportation.



TECHNICAL SPECIFICATION

Stainless Steel Grade	304 (316 also available)
Size	1200mm L x 760mm D x 900mm H
Construction	Fully – welded

FEATURES

- Strong
- High quality
- 304 and 316 grade stainless steel
- Suitable for cleanroom, laboratory and food preparation areas