

TECHNICAL DATA SHEET

Stainless Steel Wall Mounted Knee Operated Bin Bag Holder



DESCRIPTION

Wall mounted bin bag holder manufactured from grade 304 stainless steel, complete with bag clamping band and knee-operated pedal for hygienic waste disposal in cleanrooms, laboratories and other clean environments. Standard 350mm diameter. Manufactured with pre-drilled holes for easy wall mounting and stainless steel lid.



TECHNICAL SPECIFICATION

Manufacture	Optimum Protection
Material	Stainless steel
Sterility	Non-sterile
Diameter	350mm (standard size)

FEATURES

- 304 Grade Stainless Steel
- Can be manufactured in 316 grade stainless steel
- Pre-drilled holes for wall mounting
- Bag clamping band
- Knee operated pedal
- Stainless steel lid