

TECHNICAL DATA SHEET

Sterile DE Alcohol 70% - 5 Litre Bottles



DESCRIPTION

70% Denatured ethanol alcohol (DE) in water for injection quality water (WFI). All products use European Pharmacopeia grade alcohol blended with WFI quality water to guarantee endotoxin levels below 0.25 EU/ml, prior to being sterilised by a fully validated gamma irradiation process at no less than 25kGy. A premium quality product suitable for use in all cleanrooms. Double bagged and suitable for transfer.



TECHNICAL SPECIFICATION

Sterility	Sterile (Gamma Irradiated)
Filtration	0.2 micron
Size	5L Bottle
Quantity	2 x 5L bottles
Shelf Life	Unopened - 2 years from manufacture In-Use - 13 Weeks
Endotoxin Levels	<0.25 EU/ml

All products are available with full Quality Assurance documentation including:

- Certificate of Analysis – C of A
- Certificate of Irradiation – C of I
- Certificate of Conformity – C of C

FEATURES

- Fully assembled, sterile and ready for use
- Trigger sprays are individually double bagged
- Validated In-Use sterility for 13 weeks
- 0.2 micron filtered and 25kGy gamma irradiated for sterility
- Low Endotoxin level – suitable for product contact areas
- Filled and double bagged under Grade A conditions in a Grade C (ISO 7) cleanroom
- Batch coded with Detex dot to clearly indicate irradiation