

TECHNICAL DATA SHEET

Sterile DE Alcohol 70% - Trigger Spray



DESCRIPTION

70% Denatured ethanol alcohol (DE) in water for injection quality water (WFI). All products use European Pharmacopeia grade alcohol blended with WFI quality water to guarantee endotoxin levels below 0.25 EU/ml, prior to being sterilised by a fully validated gamma irradiation process at no less than 25kGy. A premium quality product suitable for use in all cleanrooms. Double bagged and suitable for transfer.



TECHNICAL SPECIFICATION

Sterility	Sterile
Product	0.2 micron filtered
Size	900ml trigger spray
Quantity	Case of 6 x 900ml bottles
Shelf Life	2 years unopened from manufacture
Compatibility	Grade A conditions in a Grade C (ISO 7) cleanroom

All products are available with full Quality Assurance documentation including:

- Certificate of Analysis – C of A
- Certificate of Irradiation – C of I
- Certificate of Conformity – C of C

FEATURES

- Fully assembled, sterile and ready for use
- Trigger sprays are individually double bagged
- Validated In-Use sterility for 13 weeks
- 0.2 micron filtered and 25kGy gamma irradiated for sterility
- Low Endotoxin level – suitable for product contact areas
- Filled and double bagged under Grade A conditions in a Grade C (ISO 7) cleanroom
- Batch coded with Detex dot to clearly indicate irradiation