

# TECHNICAL DATA SHEET

## Sterile Isopropanol (IPA) 70% 5 Litres



### DESCRIPTION

70% Isopropanol or isopropyl alcohol (IPA) in water for injection quality water (WFI). All products are manufactured in ISO Class 7 Cleanroom, and filtered to 0.2um at point of fill, double bagged and gamma irradiated 20-45kGy. A premium quality product suitable for use in all cleanrooms. Double bagged and suitable for transfer. This product is sterile and ready to use, it has a 2 year unopened shelf life. It is 0.2 micron filtered and Gamma irradiated for sterility (20-45kGy). Low endotoxin level and suitable for contact areas. Filled and double bagged under Grade A conditions in a Grade C (ISO 7) cleanroom batch coded with Detex dot to clearly indicate irradiation.



### TECHNICAL SPECIFICATION

Irradiation Standard	Gamma (20-45kGy)
Product	0.2 micron filtered
Manufacture	Agma
Size	5 Litre Bottles
Quantity	2x5 Litre bottles per case

All products are available with full Quality Assurance documentation including:

- Certification of Analysis– C of A
- Certificate of Irradiation– C of I
- Certificate of Conformity– C of C

### FEATURES

- Sterile and ready for use
- 2 x 5 Litre bottles per case
- 2 Years unopened shelf life
- 0.2 micron filtered and Gamma irradiated for sterility (20-45kGy).
- Low Endotoxin level– suitable for product contact areas
- Filled and double bagged under Grade A conditions in a Grade C (ISO 7) cleanroom
- Batch coded with Detex dot to clearly indicate irradiation