

# TECHNICAL DATA SHEET

## Sterile Latex Gloves



### DESCRIPTION

Sterile Latex Powder Free gloves for minimal risk applications. They offer superior protection while maintaining strength, sensitivity, and dexterity. The textured surface improves grip, and unlike some powder-free gloves, they are simple to put on. They have standard cuffs and are ambidextrous. Glove chlorination provides advantages such as easy glove donning and a significant reduction in the amount of extractable natural proteins present in the gloves. Absence of powder reduces the risk of powder-induced contamination. Individually wrapped in pairs to allow for easy dispensing.



### TECHNICAL SPECIFICATION

Sterility	Sterile
Quantity	Pack of 50 gloves or 1 carton – 4 packs x 50
Sizes	S (6-7) M (7-8) L (8-9) XL (9-10)
Test Standards	EN 455 Parts 1, 2, 3 & 4 AQL 1.5

### FEATURES

- Ambidextrous design
- Beaded Cuff
- Micro Textured
- Complies with EN 455, Parts 1, 2, 3 & 4
- AQL 1.5
- Smooth, Soft and Easy to Fit
- Chlorinate, Powder-Free Natural Latex
- Sterile and Bagged in Pairs
- NOTE: Not recommended for use with corrosive liquids

### CERTIFICATION

