

TECHNICAL DATA SHEET

Sterile Neutral Detergent Wipes



DESCRIPTION

Low particulate wipes impregnated with a neutral detergent in high purity water, specifically developed for cleanroom use prior to disinfection. With highly effective surface cleaning properties, it is ideal for deep cleaning prior to disinfection, to remove visible soiling or build up of residue. Manufactured using WFI quality water guaranteeing a low endotoxin level and suitability for product contact areas. 0.2 micron filtered and 25kGy gamma irradiated for sterility. Optimised fluid dosing in conjunction with the highly absorbent 68 gsm polyester/cellulose blend enables the product to hold sufficient fluid for in-use requirements without drying out too quickly.



TECHNICAL SPECIFICATION

Irradiation Standard	Gamma 25kGy
Product	0.2 micron filtered
Quantity	Pack of 10 wipes Case of 10 x 10 wipes
Wipe Dimensions	Each wipe is 28 x 38cm
Shelf Life	2 years from manufacture

FEATURES

- Low particulate wipes
- Low endotoxin level
- Highly absorbent 68 gsm polyester/cellulose wipes
- Wipes quarter folded
- Double bagged and batch coded with Detex dot
- Packed in a triple laminate, re-sealable foil bag