

TECHNICAL DATA SHEET

Sterile Nitrile Gloves



DESCRIPTION

Sterile Blue Nitrile Gloves are designed for use in clinical procedures. Highly effective and a durable barrier against blood/ blood and bodily fluids. They are perfect for use in the aid of contamination between caregivers and patients. Because the sterile gloves are individually wrapped in pairs, this allows for easy dispensing around a ward, care home or in a cleanroom environment. A convenient alternative to Latex and are comfortable and pleasant to wear. Textured throughout to provide excellent grip with a range of sizes to suit your hand shape.

Typical applications include medical, care homes, food



TECHNICAL SPECIFICATION

Sterility	Sterile
Size	S (6-7) M (7-8) L (8-9) XL (9-10)
Quantity	Pack of 50 gloves or 1 Carton– 4 packs x 50

FEATURES

- Ambidextrous design
- Beaded Cuff
- Accelerator and Powder Free
- Complies with EN 455, Parts 1, 2, 3 & 4
- AQL 1.5
- Smooth, Soft and Easy to Fit
- Sterile and Bagged in Pairs