

TECHNICAL DATA SHEET

Surgical Face Mask With Ear Loops



DESCRIPTION

Non irritating 3ply, spunbond polypropylene mask that provides high quality filtration from airborne particles and offers easy breathability. Includes ear loops. Ideal for hospital or clinical environments. Suitable for a range of procedures including medical, food processing, laboratories and cleanrooms. This mask is not classed as PPE. Fluid-resistant (Type IIR) surgical mask is Recommended for primary, outpatient, community and social care by setting, NHS and independent sector.



TECHNICAL SPECIFICATION

Manufacture	Optimum Protection
Quantity	Pack of 50 Carton of 40 x50
Material	Outer layer: Non-Woven Polypropylene Middle layer: High efficiency blown melt. Inner layer: Polyester non-woven fabric
Test Standard	EN 14683:2019 (Type IIR)
Filtration Efficiency	EN14683:2019 Annex B >98%
Breathability	<60 Pa/cm ²
Microbial Cleanliness	<30 CFU/g

FEATURES

- Size: 17cm x 9cm
- 3 Ply
- Ear loop: Urethane Elastic Fiber.
- Nose clamp: Metal wire with plastic coating
- Type IIR Fluid Repellent
- High bacterial filtration efficiency (BFE) ≥ 98%
- CE certified