

TECHNICAL DATA SHEET



Three Tier Stainless Steel Trolley – Self Assembly

DESCRIPTION

Our high quality Stainless Steel Self Assembly Trolleys are ideal for transporting equipment in cleanrooms, laboratories and other environments that require a focus on hygiene levels.



TECHNICAL SPECIFICATION

Material	Grade 201 (tube) and 430 (tray) stainless steel
Sizes	710mm (H) x 145mm (W) x 430mm (D)– Flat packed 85.5cm (L) x 53.5cm (D) x 93.0cm (H)
Weight	11kg
Sterility	Non-sterile
Load Capacity	150kg

FEATURES

- Stainless Steel
- 4 Swivel Castors
- 2 Lockable Castors
- 3 Shelves
- Simple Self Assembly
- Load Capacity 150 kg
- Comes with 4 quiet swivel castors
- Two lockable castors
- Trolley protecting disc bumpers
- Sound absorbing lining
- Comes with a 1 year warranty
- Available with 2 or 3 shelves