

TECHNICAL DATA SHEET

Tristel Duo Wipes



DESCRIPTION

Dry wipes for use with the Tristel Duo for Ultrasound disinfectant foam. Designed specifically for the even application of Tristel Duo foam, Duo Wipes are made of a soft, lint-free cloth, preventing fibre residue and scratching of delicate surfaces. Due to their low absorption power, Duo Wipes enable the optimum delivery of disinfectant onto the surface. Manufactured from 100% polypropylene, these wipes are non-woven and low-linting, and they are pre-cut to prevent shedding when dispensed.



TECHNICAL SPECIFICATION

Material	100% Polypropylene non-woven and low-linting
Thickness	18gsm
Quantity	Tubs of 200 Case of 6 x 200

STANDARDS AND DIRECTIVES:

Duo Wipes are a Class I medical device in accordance with the EU Medical Devices Regulation (MDR 2017/745)

FEATURES

- For use with Tristel Duo for Ultrasound
- 100% polypropylene
- Non-woven
- Low-linting
- Pre-cut
- 18gsm Thickness
- Minimum contact time of 30 seconds
- Tubs of 200