

# TECHNICAL DATA SHEET

## Two Tier Medical Instrument Trolley



### DESCRIPTION

Our high quality, robust, Stainless Steel medical instrument trolleys are fully welded to remove the need for self-assembly. They are ideal for use in hospital stores and operating theatres, dental laboratories and cleanroom or laboratory use. Trolley is DP1 finish as standard (brush finish). Other options include 316 Grade Stainless Steel, mirror finish and Electropolished Finish. We can manufacture any size and we can also offer instrument storage.



### TECHNICAL SPECIFICATION

|               |   |
|---------------|---|
| Material      | 304 grade stainless steel   |
| Sizes         | 45cm W x 45cm D x 87cm H<br>60cm W x 60cm D x 87cm H  |
| Load Capacity | 128kg   |
| Construction  | Fully Welded  |
| Other Options | <ul style="list-style-type: none"><li>•316 Grade Stainless Steel</li><li>•Mirror Finish</li><li>•Electropolished Finish</li><li>•Bespoke Dimensions</li><li>•Instrument Storage</li></ul> |

### FEATURES

- 304 Grade Stainless Steel - DP1 finish as standard
- 870mm to top shelf level + 20mm upturned edge
- 4 Swivel Castors (2 lockable)
- 2 Shelves
- Fully Welded Construction
- Upturned shelf edges on 3 sides. Front edge downturned
- Load Capacity : 128 kg