

TECHNICAL DATA SHEET



Two Tier Stainless Steel Trolley – Self Assembly

DESCRIPTION

Our high quality Stainless Steel Self Assembly Trolleys are ideal for transporting equipment in cleanrooms, laboratories and other environments that require a focus on hygiene levels. They require self assembly.



TECHNICAL SPECIFICATION

Material	Grade 201 (tube) and 430 (tray) stainless steel
Sizes	710mm (H) x 145mm (W) x 430mm (D)– Flat packed 84cm L x 50.5cm D x 94.5cm H
Weight	11kg
Sterility	Non-sterile
Load Capacity	150kg

FEATURES

- Stainless Steel
- 4 Swivel Castors
- 2 Lockable Castors
- 2 Shelves
- Simple Self Assembly
- Load capacity 150kg
- Comes with 4 quiet swivel castors
- Two lockable castors
- Trolley protecting disc bumpers
- Sound absorbing lining
- Comes with a 1 year warranty