

TECHNICAL DATA SHEET

Two Tier Stainless Steel Trolley – Welded



DESCRIPTION

Our high quality Stainless Steel trolleys are fully welded to remove the need for self-assembly. They are ideal for transporting equipment in cleanrooms, laboratories and other environments that require a focus on hygiene levels.



TECHNICAL SPECIFICATION

Sterility	Non-sterile
Size	90cm(L) x 60cm(W) x 95cm(H)
Load Capacity	128kg

FEATURES

- Stainless Steel
- 4 Swivel Castors
- 2 Lockable Castors
- 2 Shelves
- Fully Welded Construction
- Load Capacity : 128 kg