



# PRODUCT DATA SHEET

Product Code: **BY750S**

## Byotrol Trionic Biocide Spray

- Alcohol-Free
- Quat-Based Biocide
- Cleans and Disinfects
- Ready to Use (RTU)
- Also Available in Pre-Wetted Wipes

A non-alcohol formulation for the decontamination of medical devices and equipment, cleaning surfaces (including the removal of bio-film), and then disinfecting with broad spectrum biocidal efficacy.

The ready-to-use spray is ideal for quick cleaning and disinfection of the surfaces of medical devices in clinical and other human health environments, as well as for use in cleanrooms and laboratories.

The disinfectant kills bacteria, fungi, yeasts, viruses and mycobacteria, and has been tested to be effective under a number of microbiological test conditions. For a full list, see the next page of this data sheet.

Also available in wipe format.



V1.0 07/12/2022



# PRODUCT DATA SHEET

Product Code: **BY750S**

## Microbiological Testing

Byotrol Trionic has been tested to be effective under the following microbiological test conditions:

- EN1040 (Bactericidal)
- EN1276 (Bactericidal)
- EN13727 (Bactericidal)
- EN1650 (Fungicidal)
- EN13697 (Bactericidal, Yeasticidal, Fungicidal)
- EN13624 (Fungicidal, Yeasticidal)
- EN14476 (Virucidal)
- EN14348 (Mycobacterial)

## Health & Safety

Use biocides safely. Always read the label and product information before use. Safety Datasheets are available to download from our website.

## Applications

Ensure surfaces are visibly clean.

Wipe the contaminated surface and dispose of the wipe.

Allow surfaces to dry.

For professional use only. Dispose of wipe with general waste.

V1.0 07/12/2022