



PRODUCT DATA SHEET

Product Code: MJZK19

Chemical Spill Kit

- Each Pack Absorbs 1L of Chemical Spill
- Safe
- Efficient
- Conforms to HSE & BSI Guidelines
- Ready to Use
- Easy to Keep
- 3x Single Use Packs

Multi-use spill kit for protection against hazardous chemicals and contaminating fluids. One kit contains 3 single-use packs, and each pack absorbs up to 1 litre of chemical spill. Provides safe, efficient cleaning of chemical spillages in laboratory environments.

The absorbent pad system is tested and approved to be used with over 60 different chemical spills. Each pack also contains all PPE to ensure fast user safety in the event of a hazardous chemical spill. The kit also include a wall bracket to ensure accessibility in areas where it may be needed.

The sorbent is made up of two layers that are sonically welded together. The top layer is made of a no lint, durable spun-bond polypropylene while the core/base layer is made from high loft melt-blown fibres. In addition the core/base layer has low slip treatment for areas of high foot traffic.

Conforms to HSE and BSI guidelines for cleanup of chemical spills. The small pack makes it easy to keep ready for use in all areas and ensure accessibility where it may be needed.

The one pack use format prevents fragmentation of components and ensures all items are present when opened.



V1.0 09.06.2021



PRODUCT DATA SHEET

Product Code: MJZK19

Chemical Spill Kit

Designed and tested for clean-up of the following chemicals and fluids:

- Acetaldehyde
- Acetic Acid
- Acetic Anhydride
- Acetone
- Acetyl Chloride
- Acrolein
- Acrylic Acid
- Acrylic Emulsions
- Acrylonitrile
- Allyl Alcohol
- Aminobenzoic Acid
- Ammonia (Anhydrous)
- Ammonium Hydroxide
- Amyl Acetate
- Amyl Alcohol
- Aniline
- Aqua Regia
- Aviation Fuel
- Benzene
- Benzoic Acid
- Benzointrile
- Benzoyl Chloride
- Benzyl Alcohol
- Boric Acid
- Brake Fluid
- Bromine
- Butyl Acetate
- Butyl Alcohol
- Butylamine
- Butyric Acid
- Calcium Hydroxide
- Carbolic Acid
- Carbon Disulfate
- Carbon Tetrachloride
- Castor Oil
- Chlorine Water
- Chloroacetic Acid
- Chlorobenzene
- Chloroform
- Chlorosulfonic Acid
- Chromic Acid (50%)
- Citric Acid
- Clorox (Full Strength Bleach)
- Corn Oil
- Cottonseed Oil
- Cresol
- Cyclohexane
- Detergents
- Dichlorobenzene
- Diethyl Ether
- Diethylamine
- Dinitrobenzene
- Dioxan
- Dioctyl Phthalate
- Ether
- Ethyl Acetate
- Ethyl Alcohol
- Ethyl Benzene
- Ethyl Chloride
- Ethyl Ether
- Ethyl Propionate
- Ethylene Glycol
- Formaldehyde
- Formic Acid
- Fuel Oil
- Gasoline
- Gearbox Oil
- Glacial Acetic Acid
- Glycerol
- Heptane
- Hexane
- Hydrazine
- Hydrochloric Acid
- Hydrogen Cyanide
- Hydrogen Peroxide
- Isobutyl Alcohol
- Isobutyric Acid
- Isopropyl Acetate
- Isopropyl Alcohol
- Kerosene
- Keytones
- Linseed Oil
- Lubricating Oil
- Magnesium Hydroxide
- Methyl Alcohol
- Methyl Chloride
- Methyl Ether
- Methyl Ethyl Ketone
- Methyl Propionate
- Mineral Oil
- Motor Oil
- Naphthalene
- Nitric Acid
- Nitrobenzene
- Nitrobenzoic Acid
- Nitrotoluene
- Octane
- Oleic Acid
- Olive Oil
- Paraffin
- Perchloroethylene
- Petroleum Ether
- Phenol
- Phosphoric Acid
- Plating Solutions
- Potassium Hydroxide
- Propano
- Propionic Acid
- Propyl Alcohol
- Propylene Glycol
- Quinoline
- Resorcinol
- Salt Solutions
- Silicone Oil
- Silver Nitrate
- Soap Solution
- Sodium Bicarbonate
- Sodium Chloride
- Sodium Hydroxide
- Sodium Hypochlorite
- Sodium Nitrate
- Stannic Chloride
- Starch
- Styrene
- Sucrose
- Sulfuric Acid
- Synthetic Motor Oil
- Tannic Acid
- Toluene
- Transformer Oil
- Trichloroethylene
- Triethylene Glycol
- Turpentine
- Vinegar
- Vinyl Acetate
- Xylene

V1.0 09.06.2021



PRODUCT DATA SHEET

Product Code: MJZK19

Chemical Spill Kit

Each Pack Contains:

- 1x Pair of Nitrile Gloves
- 1x Apron or Thumb Loop Gown
- 1x Face Mask / Eye Shield
- 3x IPA Surface Wipes
- 1x Disposable Bag
- 2x Absorbent Pads
- 1x Set of Tweezers
- 1x Clinical Waste Bag
- Easy to Follow Illustrated Instructions

Each Kit Contains:

- 3x MJZ019 Packs*
- 1x Durable Red Case
- 1x Wall Bracket

*MJZ019 packs are sold on our website and can be used as refills for this kit.

V1.0 09.06.2021