



PRODUCT DATA SHEET

Product Code: CWS07

Clear Domestic Waste Bags

- Clear Domestic Waste Bags
- Wide Range of Applications
- CHSA Approved
- Printed for NHS or Unprinted
- Supplied in Rolls

CHSA approved clear / transparent polythene waste bag for use in acute and non-acute healthcare settings including hospitals, clinics, GP surgeries, care homes, dental and veterinary practices where professional quality general waste bags are required. Available in a range of sizes and strengths, and with or without printing.

Suitable for general hospital and healthcare waste, household waste and domestic waste.

No suitable for use with medical or clinical waste, infectious waste, anatomical waste, placental waste, cytotoxic or cytostatic waste. Recyclable items should not be placed in general waste if an appropriate recycling stream is available.

Bags Available:

Small Printed - 535 x 595mm, 15mu, 30L

Large Printed - 711 x 960mm, 15mu, 90L

Unprinted Light Duty - 736 x 965mm, 90L

Unprinted Medium Duty - 736 x 965mm, 90L, 10Kg Max. Weight (CHSA)

Unprinted Heavy Duty - 736 x 965mm, 90L, 15Kg Max. Weight (CHSA)



V1.0 09.06.2021