



PRODUCT DATA SHEET

Product Code: CHEM-W-I

Instrument Disinfectant Wipes

- Powerful Disinfectant
- Tubs or Pouches
- Nonwoven Wipes
- Alcohol Free
- Effective Against a Range of Micro-Organisms
- Non-Corrosive
- Class IIb CE Marked

High level non-woven disinfectant wipes for the cleaning of non-surgically invasive instruments, such as transvaginal and transrectal probes. Also suitable for use on hard surfaces of non-invasive, non-lumened medical devices, such as monitors, radiographic and diagnostic equipment.

The wipes are impregnated with a multi-active disinfectant solution that is alcohol free and tested to EN standards as effective against a range of bacteria, viruses, fungi, yeasts and mycobacteria.

CE marked as a class IIb medical device, these wipes are developed and tested by qualified chemists and microbiologists. They are then subjected to stringent EN testing procedures to ensure efficacy against a wide range of bacteria, viruses, fungi and spores.

Not to be used to disinfect surgically invasive medical devices.



V1.0 02.06.2021



PRODUCT DATA SHEET

Product Code: CHEM-W-I

This product has been tested effective against:

Virucidal Activity

- Hepatitis B
- Herpes Simplex
- Norovirus
- Vaccinia
- Adenovirus

Bactericidal Activity

- *Pseudomonas aeruginosa*
- *Escherichia coli*
- *Staphylococcus aureus*
- *Enterococcus hirae*
- *Legionella pneumophila*
- *Salmonella typhimurium*
- *Listeria monocytogenes*
- Methicillin-resistant *Staphylococcus aureus*
- Methicillin-resistant *Staphylococcus pseudointermedius*

Bacteriophage Activity

- *Lactococcus lacti*

Fungicidal Activity

- *Candida albicans*
- *Trichophyton rubrum*

Compliance

This product has been tested against the following standards:

- EN 1276
- EN 13697
- EN 1650
- EN 14476

V1.0 02.06.2021