

PRODUCT SPECIFICATION

Product Name :	Isopropanol			
Specification Reference :	IS003\11			
Issue No :	2			
Date of Issue :	23/11/9			
Properties	Units	Minimum	Maximum	Test Method
Purity	%(w/w)	99.8		DIN 55685
Water	%(w/w)		0.1	ASTM D1364
Acidity as acetic acid	%(w/w)		0.001	ASTM D1613
Water Miscibility		Miscible		ASTM D1722
Non-volatile matter	g/100 ml		0.001	ASTM D1353
Appearance	Clear and free from suspended matter			ASTM D4176
Colour, Pt-Co			5	ASTM D1209 (4)
Density @ 20°C	g/ml	0.785	0.786	ASTM D4052 (4)
Distillation IBP	°C	81.8		ASTM D1078 (4)
Distillation DP	°C		82.8	ASTM D1078 (4)
Refractive Index @ 20°C		1.376	1.378	ASTM D1218 (4)

Additional Information:

(4) – Agreed specification limits, however, statistical average value reported on COA

Product as produced complies with ASTM D770, DIN 53245 & ACS 8th edition

This certificate has been computer generated and therefore is not signed.

Any suggestion or representation concerning any possible use of the product made or supplied by, whether in literature or in response to a specific enquiry, is given in good faith, but it is entirely for the end user to satisfy itself fully as to the suitability of the product for any particular purpose. No warranty or representation relating to such possible use is expressly given or shall be implied (whether by statute or otherwise) and shall be excluded, as far as is legally permissible.

