

## Product Description

Our high quality Stainless Steel Trolleys are fully welded to remove the need for self-assembly. They are ideal for transporting equipment in cleanrooms, laboratories and other environments that require a focus on hygiene levels




**BESPOKE**

Do you have a specific furniture design that you cannot find

Utilise our bespoke manufacturing service by contacting us using the details below

**FURNITURE**



## Key Features

- Manufactured from 304 grade Stainless Steel (See Below)
- Carrying capacity 128kg
- Product comes with a 1 year warranty
- Comes with four quiet swivel castors
- Two Lockable Castors
- Available with 2 or 3 Shelves
- Sound Insulation
- Trolley protecting disc bumpers

## Delivery Dimensions

- Product is delivered fully welded, removing the need for assembly
- Dimensions: 1060(H) x 920 (L) x 620 (D) mm
- Weight: 24kg
- Fully Welded Construction

Shelves	Length (L) (cm)	Height (H) (cm)	Depth (D) (cm)	Distance Between Shelves (cm)	
2	90	90	60	58.7	
3	90	90	60	A = 27.5	B = 27.5

### 304 Stainless Steel

304 grade is the most common form of stainless steel. It contains between 16 and 24 percent chromium and up to 35 percent nickel - as well as small amounts of carbon and manganese. 304 stainless steel has a high resistance to rust. It withstands corrosion from most oxidizing acids. However, it is susceptible to corrosion from chloride solutions (notably saline environments with high amounts of sodium chloride).

**Note:** We can also provide trolleys in 316 Stainless Steel. This would be a bespoke item. Please contact us to discuss your requirements.