

# PRODUCT DATA SHEET

## Sterile AZO Wipes

### Benefits:

- Broad spectrum activity against a range of the most 14 pathogenic micro-organisms including *E. Coli*, *Pseudomonas aeruginosa*, *Proteus mirabilis*, *Staphylococcus aureus*, *Enterococci*, *Mycobacterium terrae*, MRSA, Hepatitis B, Hepatitis C, HIV, Vaccinia, Flaviviridae, *Aspergillus niger* and *Candida albicans*
- Extremely fast contact times from **15 seconds**
- Wipes are low-linting, strong and durable for everyday use
- Wipes are non-smearing. Does not streak hard surfaces
- Sprays comes with lock mechanism to prevent leaking and fire risk
- Certified Halal approved by European Halal Development agency. Suitable for food manufacturing, hospitality & catering sectors
- Compliant with HTM-01-05 guidelines



**AZOWIPE™** has been the leading product of its type for over 30 years. **AZOWIPE™** is convenient for use in all clinical environments. With the increase in virulent organisms the introduction of **AZOWIPE™** as part of a robust infection control protocol is now even more important. Isopropyl alcohol, the active component of **AZOWIPE™** is acknowledged by leading authorities to be an effective medium for the routine disinfecting of clinical surfaces. Kills bacteria in under a minute - proven to be effective against a variety of organisms including MRSA.

AZOWIPE™ is compliant with the following recognised European EN standards against a wide range of micro-organisms including:

- Bactericidal – EN1276 and EN13727
- Virucidal – EN14476
- Yeasticidal/Fungicidal – EN1650/EN13624
- Mycobactericidal – EN14348

### Gamma Irradiation:

The treatment of the Sterile Azowipes, Gamma processing, has been one of the principal methods of sterilisation in the healthcare industry since the introduction of the concept of single use, sterile, disposable medical devices in the 1960s.

This simple, proven process is safe, reliable, and highly effective at treating single-use medical devices. With the ability to penetrate products while sealed in their final packaging ensures full sterility of the product

### Available in:

- Single Bagged
- Double Bagged
- Triple Bagged

Product	Description	Size	Pack Qty	Case Qty
AZOWIPE STERILE	70% IPA disinfectant wipes; gamma irradiated canister, sealed and sterilised in polythene bags for use in cleanrooms and other areas where aseptic technique is essential	20x22cm	200	12

